**ARGON MEDICAL DEVICES, INC.**  
1445 Flat Creek Road  
Athens, Texas 75751 USA  
Tel: 800-927-4669;  
Tel: +1 903-675-9321  
www.argonmedical.com


**Emergo Europe, Molenstraat 15, 2513 BH.**  
The Hague, The Netherlands  
Tel: +31 (0)70 345 8570  
Fax: +31 (0)70 346 7299

STERILE EO

0086

## ENGLISH

**Indications for Use:** Intended for the purpose of harvesting bone and/or bone marrow specimens.

**Contraindications:** For use only for biopsies of bone/bone marrow as determined by a licensed physician. These needles should be used by a physician familiar with the possible side effects, typical findings, limitations, indications and contraindications of bone biopsy. Physician judgment is required when considering biopsy on patients with bleeding disorder, or receiving anti-coagulant medications.

**Caution:** For Single Patient Use Only. Do not attempt to clean or resterilize this product. After use, this product may be a potential biohazard. Handle in a manner which will prevent accidental puncture. Dispose in accordance with applicable laws and regulations.

**Note:** These instructions are NOT meant to define or suggest any medical or surgical technique. The individual practitioner is responsible for the proper procedure and techniques to be used with this device.

### Aspiration Procedure:

- Place the patient in a right or left lateral position, with the back comfortably flexed and the top knee drawn toward the chest.
- Locate the posterior superior iliac spine and mark.
- Using sterile technique, prepare the skin with antiseptic and drape.
- Infiltrate the marked area with local anesthetic, especially the periosteum.
- Make a skin incision with a scalpel blade over the marked area.
- Hold the needle with the proximal end in palm and the index finger against the shaft near the tip. This position stabilizes the needle and allows better control.
- Introduce the needle through the incision pointing toward the anterior superior iliac spine and bring it into contact with the posterior iliac spine.
- Using gentle, but firm pressure, advance the needle tip through the periosteum, into the cortex, by rotating the needle in an alternating clockwise-counterclockwise motion.
- Remove stylet by rotating upper section of handle 90°, and pulling straight out.
- Continue advancing the needle cannula forward, while rotating in an alternating clockwise-counterclockwise motion, into the bone marrow cavity. Entrance into the marrow cavity is generally detected by decreased resistance.
- Attach a syringe with a luer taper to the hub of the Bone Marrow Biopsy/Aspiration Needle with a firm, push, twist motion.
- Apply suction by withdrawing syringe plunger. Remove the aspirated specimen collected from the syringe.

### Biopsy Procedure:

Take the biopsy after the aspiration using the same skin incision, but choosing a different location to obtain the biopsy sample.

Repeat Steps 1 through 10 of the procedure and continue as follows:

- Optional - If desired, cannula can be capped with luer cap provided.
- Slowly and gently advance the needle a millimeter at a time, with clockwise-counterclockwise motion (for better cutting) until adequate marrow is obtained.
- To harvest the specimen:
  - With negative pressure,
    - Insert the probe into the cutting cannula to check the sample length in the needle lumen.
    - Attach a syringe to the fitting on the hub on the cannula handle and draw negative pressure to hold your specimen.
    - Remove needle.
    - Push the specimen out with the probe.
  - Mechanically with the use of the T-Lok™ Extraction cannula,
    - Insert the probe into the cutting cannula to check the sample length in the needle lumen.
    - Remove the probe and insert the T-Lok™ Extraction cannula FULLY into the needle cannula.
    - Slowly rotate and remove the needle and T-Lok™ Extraction cannula TOGETHER.
    - Remove the T-Lok™ Extraction cannula from the needle cannula and push the specimen out with the probe.

**Note:** For safety and convenience, Medical Device Technologies has introduced the Anti-Stick Probe Guide as a standard feature of the "J" type Bone Marrow Needle. The purpose of this device is two-fold: provide an easier method to align the probe in the needle tip for sample expulsion and protect the sharp tip to eliminate any accidental injury during sample expulsion.

After sample is obtained, slip larger opening end of guide over the needle tip. Push guide down on needle until it fits snugly. Insert probe into smaller opening and push gently to expel sample from proximal end of needle. Remove probe from guide. If another sample is to be taken, remove guide and reassemble needle stylet and cannula. Wipe cannula and stylet tip clean. If no other samples are to be obtained, leave guide in place to cover tip of needle and discard needle assembly.

T-Lok™ is a trademark of Medical Devices Technologies, Inc.

**Eng: T-Lok™ Bone Marrow Biopsy Needle**  
**Ger: T-Lok™ Knochenmark-Biopsienadel**  
**Fre: Pour biopsie de la moëlle osseuse T-Lok™**  
**Spa: Aguja de biopsia para médula ósea T-Lok™**  
**Ita: Ago per biopsia del midollo osseo T-Lok™**  
**Swe: T-Lok™ Biopsinål för benmärg**  
**Por: Agulha para Biópsia de Medula Óssea T-Lok™**  
**Gre: Βελόνα βιοψίας Μελού Οστού T-Lok™**  
**Dut: T-Lok™ beenmergbiopsienaald**  
**Fin: T-Lok™ luuydinbiopsineula**  
**Nor: T-Lok™ Biopsinål for benmarg**  
**Dan: T-Lok™ Kanyle til Knoglemarvsbiopsi**  
**Tur: T-Lok™ Kemik iliği Biyopsisi İğnesi**

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